

Training Manual for PISTA Alzheimer's Table Tennis

Introduction

Alzheimer's disease is a chronic and progressive brain disorder that primarily affects individuals over the age of 65, although early-onset Alzheimer's can occur in younger adults. This devastating condition is characterized by the accumulation of abnormal protein deposits, such as beta-amyloid plaques and tau tangles, in the brain. These accumulations interfere with the normal functioning of neurons, leading to their damage and death.

The hallmark of Alzheimer's is a decline in cognitive function, with memory loss being one of the most prominent early symptoms. As the disease advances, individuals may struggle with language, have difficulty recognizing familiar faces, experience mood swings and personality changes, and ultimately lose the ability to perform even the most basic daily tasks like dressing, eating, and bathing.

Alzheimer's places an immense burden not only on those directly affected but also on their families and caregivers who often provide round-the-clock support. As a result, the disease has significant emotional, financial, and societal implications. While there is currently no cure for Alzheimer's, ongoing research aims to better understand its underlying mechanisms and develop more effective treatments to alleviate its symptoms and potentially slow down its progression.

In summary, Alzheimer's disease is a complex and deeply impactful neurological condition that robs individuals of their cognitive abilities, memory, and ultimately, their independence. It poses a substantial challenge to healthcare systems and families worldwide, making it a subject of critical importance in the field of medical research and public health.

There is currently no cure for AD, however, there are treatments that can slow the progression of the disease and improve quality of life for patients. Education on this disease is important for physiotherapists because in Canada, there are an estimated 747,000 people diagnosed with AD and Dementia per year. As the population of older adults increase, this number is projected to reach 1.4 million in 2031.

Physiotherapists play an important role in treating patients with AD since exercise has been shown to improve physical function, functional independence, and slow the progression of the disease. This case study outlines the role that physiotherapy plays in the assessment, treatment, and follow-up of a patient with AD.

One of the challenges of dealing with the patient in this case is recognizing that his memory loss and the progressive nature of the disease will often lead to worsening of symptoms over time. The physiotherapist took steps to ensure the patient's understanding of treatment

measures and kept in communication with his wife to ensure proper delivery of the intervention at home.

Behavioral and Psychological Symptoms of Dementia (BPSD)

Behavioural and Psychological Symptoms of Dementia (BPSD) will develop in more than 90% of individuals diagnosed with dementia. Symptoms may include delusions, hallucinations, aggression, screaming, restlessness, wandering, depression, and anxiety. Other symptoms include disinhibition, sexual behaviours, apathy, sleep disturbances, and compulsive or repetitive behaviours. BPSD results in impaired quality of life, increased cost of care, rapid cognitive decline, and massive caregiver burden

Behavioral and psychological symptoms of dementia (BPSD), also known as neuropsychiatric symptoms, represent a heterogeneous group of non-cognitive symptoms and behaviors occurring in subjects with dementia. BPSD constitute a major component of the dementia syndrome irrespective of its subtype. They are as clinically relevant as cognitive symptoms as they strongly correlate with the degree of functional and cognitive impairment. BPSD include agitation, aberrant motor behavior, anxiety, elation, irritability, depression, apathy, disinhibition, delusions, hallucinations, and sleep or appetite changes. It is estimated that BPSD affect up to 90% of all dementia subjects over the course of their illness, and is independently associated with poor outcomes, including distress among patients and caregivers, long-term hospitalization, misuse of medication, and increased health care costs. Although these symptoms can be present individually it is more common that various psychopathological features co-occur simultaneously in the same patient.

BPSD symptoms can be organized into 5 symptom clusters:

1. **Apathy** (lack of initiative)
2. **Psychosis** (delusions, hallucinations)
3. **Aggression/Agitation** (verbal, physical)
4. **Hyperactivity** (pacing, restlessness, disinhibition)
5. **Affective** (dysphoria, elation irritability, anxiety)

How does Alzheimer's disease affect the brain?

The ongoing research into Alzheimer's disease has been uncovering the intricate and multifaceted changes that occur within the brain. Remarkably, these changes can commence a decade or more before any noticeable symptoms manifest. In this early, preclinical stage of Alzheimer's, the brain undergoes a series of distressing transformations, marked by the accumulation of abnormal proteins, particularly amyloid plaques and tau tangles. These protein aggregates interfere with the normal functioning of neurons, causing previously healthy nerve cells to malfunction, lose their connections with neighboring cells, and eventually perish.

While the formation of amyloid plaques and tau tangles is a prominent feature of Alzheimer's pathology, it represents just one facet of a much larger and intricate puzzle. Numerous other complex brain changes are believed to be involved in the development and progression of the disease, further deepening our understanding of its complexity.

Initially, the damage appears to concentrate in crucial brain regions, notably the hippocampus and the entorhinal cortex, both vital for memory formation and storage. As the disease advances, the destructive process spreads beyond these regions, affecting additional parts of the brain. This widespread impact leads to observable brain atrophy, or shrinkage, as neurons continue to degenerate and die.

By the final stage of Alzheimer's, the cumulative damage is extensive, and brain tissue has undergone significant reduction. This widespread devastation underscores the severity of the disease and its profound impact on an individual's cognitive abilities, memory, and overall brain function. As researchers delve deeper into these intricate brain changes, they hold hope for uncovering novel therapeutic strategies to address the underlying mechanisms of Alzheimer's and, ultimately, to find more effective treatments for this devastating condition.

Signs of Alzheimer's Disease

Alzheimer's symptoms can vary greatly among individuals. Usually, memory problems are one of the signs of the disease. However, in the stages of Alzheimer's, other aspects of thinking and understanding may also be affected. This includes having difficulty finding the words, understanding images and spatial relationships and making judgments. As Alzheimer's progresses these symptoms become more pronounced. Can lead to increased confusion and changes, in behavior.

For people with Alzheimer's those with late onset Alzheimer's symptoms typically appear in their mid-60s or later. In cases where the disease begins before age 65 it is known as early onset Alzheimer's and can start early as a person's 30s.

From a clinical perspective Alzheimer's usually advances through different stages; preclinical stage (before noticeable symptoms) mild stage (also called early stage) moderate stage and severe stage (also referred to as late stage). These stages help healthcare professionals and caregivers understand and address the evolving challenges faced by individuals, with Alzheimer's and their families as the disease progresses.

Preclinical Alzheimer's disease

Scientific research indicates that the intricate alterations in the brain linked to Alzheimer's disease, such as the formation of amyloid plaques and tau tangles, commence approximately ten or more years before observable memory and cognitive difficulties emerge. This preliminary phase, characterized by these brain changes preceding the onset of dementia, is referred to as preclinical Alzheimer's. Nevertheless, it is essential to emphasize that not all individuals experiencing these brain changes ultimately progress to dementia.

Signs of Mild Alzheimer's disease

In mild Alzheimer's, a person may seem healthy but has more and more trouble making sense of the world around them. The realization that something is wrong often comes gradually to the person and their family. Problems can include:

- Memory loss that disrupts daily life
- Poor judgment, leading to bad decisions
- Loss of spontaneity and sense of initiative
- Losing track of dates or knowing current location
- Taking longer to complete normal daily tasks
- Repeating questions or forgetting recently learned information
- Trouble handling money and paying bills
- Challenges in planning or solving problems
- Wandering and getting lost
- Losing things or misplacing them in odd places
- Difficulty completing tasks such as bathing
- Mood and personality changes
- Increased anxiety and/or aggression

Alzheimer's is often diagnosed at this stage.

Signs of moderate Alzheimer's disease

In this stage, more intensive supervision and care become necessary. These changes and increasing needs can be difficult for many spouses and families. Symptoms may include:

- Increased confusion and memory loss, such as forgetting events or personal history
- Withdrawal from social activities
- Inability to learn new things
- Difficulty with language and problems with reading, writing, and working with numbers
- Difficulty organizing thoughts and thinking logically
- Shortened attention span
- Problems coping with new situations
- Changes in sleeping patterns, such as sleeping more during the day and being restless at night
- Difficulty carrying out familiar, multistep tasks, such as getting dressed
- Occasional problems recognizing family and friends
- Hallucinations, delusions, and paranoia
- Impulsive behavior, such as undressing at inappropriate times or places, or using vulgar language
- Inappropriate emotional outbursts
- Restlessness, agitation, anxiety, tearfulness, wandering — especially in the late afternoon or evening
- Repetitive statements or movement, occasional muscle twitches

Signs of severe Alzheimer's disease

People with severe Alzheimer's cannot communicate and are completely dependent on others for their care. Near the end of life, the person may be in bed most or all of the time as their body shuts down. Symptoms often include:

- Inability to communicate
- No awareness of recent experiences or surroundings
- Weight loss with little interest in eating
- Seizures

- General physical decline, including dental, skin, and foot problems
- Difficulty swallowing
- Groaning, moaning, or grunting
- Increased sleeping
- Loss of bowel and bladder control

A common cause of death for people with Alzheimer's is aspiration pneumonia. This type of pneumonia develops when a person cannot swallow properly and takes food or liquids into the lungs instead of air.

Symptoms of mild cognitive impairment

Some people have a condition called mild cognitive impairment (MCI), which can be an early sign of Alzheimer's. However, not everyone with MCI will develop Alzheimer's. People with MCI can still take care of themselves and perform their normal activities. MCI memory problems may include:

- Losing things often
- Forgetting to go to events or appointments
- Problems communicating because of difficulty finding words

Main Characteristics of The Alzheimer's Disease

Alzheimer's disease is characterized by a range of cognitive and behavioral symptoms that progressively worsen over time. The main characteristics of Alzheimer's include:

1. **Memory Loss:** One of the earliest and most prominent signs of Alzheimer's is the loss of short-term memory. Individuals may have difficulty remembering recent events, appointments, or conversations. As the disease advances, long-term memories may also become affected.
2. **Cognitive Decline:** Alzheimer's leads to a decline in various cognitive functions, including problem-solving, decision-making, attention, and reasoning. Individuals may struggle with tasks that were once routine.

3. **Language Problems:** People with Alzheimer's often have difficulty finding the right words, following or joining in conversations, and may repeat themselves frequently. Writing and reading skills can also deteriorate.
4. **Disorientation:** Individuals may become disoriented regarding time and place, leading to confusion about their current location or the date. This can cause anxiety and agitation.
5. **Impaired Judgment:** Alzheimer's can impair a person's judgment and decision-making abilities, potentially leading to poor financial choices or unsafe behavior.
6. **Personality and Behavior Changes:** Behavioral changes are common, including mood swings, irritability, aggression, depression, and apathy. Individuals may also exhibit socially inappropriate behavior.
7. **Difficulty with Routine Tasks:** As the disease progresses, individuals may struggle with everyday tasks like dressing, eating, and bathing, eventually requiring assistance with these activities.
8. **Loss of Independence:** Alzheimer's disease gradually erodes a person's independence, making them increasingly reliant on others for care and support.
9. **Changes in Spatial Abilities:** Individuals may have trouble with spatial tasks, such as getting lost in familiar surroundings or being unable to recognize familiar faces.
10. **Inability to Recognize Loved Ones:** In advanced stages, people with Alzheimer's may struggle to recognize family members and close friends.
11. **Motor and Coordination Problems:** Some individuals may experience motor difficulties, such as problems with coordination and balance.
12. **Wandering and Agitation:** Restlessness, pacing, and wandering are common behavioral symptoms, which can be challenging for caregivers to manage.
13. **Sleep Disturbances:** Sleep problems, including insomnia and daytime drowsiness, are often seen in Alzheimer's patients.

It's important to note that Alzheimer's disease progresses differently in each individual, and the severity and timeline of these symptoms can vary widely. Early diagnosis and intervention are crucial to manage the condition and provide the best possible quality of life for affected individuals.

Mechanisms of Memory Loss in Alzheimer's Disease

Memory loss in Alzheimer's disease is a central and often deeply distressing aspect of the condition. It represents a gradual and relentless decline in an individual's ability to remember,

affecting various forms of memory, from short-term to long-term, and encompassing a range of cognitive processes. At the core of this memory impairment lies the progressive deterioration of brain regions critical for memory, notably the hippocampus and its surrounding structures.

Early on, individuals with Alzheimer's may exhibit subtle signs of forgetfulness, such as misplacing everyday objects, forgetting appointments, or repeating questions. These initial lapses in short-term memory can be attributed to the hippocampus's diminished capacity to encode and retrieve new information effectively. As the disease advances, this difficulty in forming new memories becomes more pronounced, making it challenging to remember recent events, conversations, or even the names of close family members.

Long-term memories, stored in regions of the brain outside the hippocampus, tend to remain relatively preserved in the early stages. However, as Alzheimer's inexorably spreads throughout the brain, even these cherished memories may begin to fade. Individuals may struggle to recall significant life events, important personal details, or historical facts that they once knew effortlessly.

The impact of Alzheimer's-related memory loss is not limited to cognitive functioning alone. It profoundly disrupts daily life, diminishing one's ability to perform routine tasks, maintain personal relationships, and engage in meaningful activities. It also introduces a considerable emotional burden, as individuals with Alzheimer's often become aware of their declining memory, leading to frustration, anxiety, and even depression.

Family members and caregivers also face significant challenges in coping with the progressive memory loss of their loved ones. They may need to provide increasing levels of support and supervision as their relative's ability to manage daily activities independently declines.

It is important to emphasize that while Alzheimer's-related memory loss is devastating, early diagnosis and timely interventions, such as medication and cognitive therapies, can help manage the symptoms and enhance the quality of life for affected individuals. Moreover, ongoing research into Alzheimer's disease aims to develop treatments that can potentially slow or halt the progression of memory loss and other cognitive symptoms, offering hope for a brighter future in the fight against this debilitating condition.

Brain waves clear Alzheimer's plaques

Our brains can get quite rhythmical as groups of neurons fire together in oscillating patterns during various activities, such as processing sensory information or encoding memories. A study suggests that one of these rhythmical patterns, called gamma waves, can stimulate the brain's immune cells to clear out amyloid- β ($A\beta$), the peptide that forms plaques in the brains of people with Alzheimer's disease (Nature 2016, DOI: 10.1038/nature20587).

The researchers, led by Li-Huei Tsai of Massachusetts Institute of Technology, think their results hint at a noninvasive therapy for the neurodegenerative disorder.

"The implications are significant," says Michal Schwartz of the Weizmann Institute of Science, who was not involved in the work. And, she says, if the effects of gamma waves also improve cognition in Alzheimer's models, "it's unbelievable."

Gamma waves—patterns of brain activity with frequencies between 30 and 90 Hz—are often disrupted in many neurological disorders, including Alzheimer's. The MIT team wanted to understand the connection between these oscillations and the disease.

The scientists used optogenetics to generate gamma waves at 40 Hz in the brains of mice engineered to overproduce $A\beta$ and found that the rhythmical activity significantly reduced levels of the peptide. When the scientists studied these animals' brains further, they found evidence that immune cells called microglia had cleared the peptide.

But optogenetics, which activates neurons with light beamed into the brain, is invasive. To demonstrate a noninvasive way to trigger gamma waves, the team placed the animals in a dark chamber and then flashed a light-emitting diode at 40 Hz. This flickering generated gamma waves in the visual cortex of the animals' brains and reduced the amount of amyloid plaques in that brain region after a week of treatment.

Noninvasive brain wave treatment reduces Alzheimer's pathology, improves memory in mice

Recent research supported by the NIA has demonstrated that stimulating gamma brain waves in mice can reduce Alzheimer's-related proteins and slow down the neurodegeneration associated with the disease. Healthy brains exhibit various rhythmic patterns, known as brain waves, which operate at different frequencies. Gamma brain waves, which oscillate between 30 to 100 Hz, are linked to higher-order cognitive functions and tend to decrease in people with Alzheimer's.

In previous studies at the Massachusetts Institute of Technology (MIT), researchers discovered that exposing Alzheimer's mouse models to LED lights flickering at 40 Hz could stimulate gamma waves. This not only lowered levels of beta-amyloid and tau, the proteins associated with Alzheimer's, but also enhanced the activity of microglia in clearing harmful debris.

In a study reported in the journal *Cell*, MIT researchers delved further into the effects of gamma wave stimulation, this time using sound in Alzheimer's mouse models. Clicks played at 40 Hz for an hour a day over a week reduced beta-amyloid levels in the auditory cortex and nearby hippocampus, which is responsible for learning and memory. The stimulated mice performed better on memory tasks, and there were notable changes in the activation responses of microglia, astrocytes, and blood vessels.

Subsequently, the researchers exposed mice to a combination of light and sound stimulation, expanding the effects beyond the visual and auditory cortex to the prefrontal cortex, an area vital for planning and task completion. Analysis revealed a unique clustering effect of microglia around amyloid deposits in the stimulated mice, along with reduced amyloid pathology. However, these effects were short-lived and diminished after a week.

In another study published in the journal *Neuron*, MIT researchers tested the effects of longer-term treatment by subjecting mouse models with advanced Alzheimer's disease to up to six weeks of gamma wave stimulation via visual methods. The results demonstrated that stimulation increased gamma brain waves in various brain regions, preserved neuronal and synaptic density, improved performance on memory tasks, and reduced inflammation. These findings indicate an overall neuroprotective effect, even in the later stages of neurodegeneration.

The results of this research build on prior investigations of gamma wave stimulation as a potential treatment for Alzheimer's disease, and researchers are now planning studies to assess whether these benefits could extend to humans.

Does Synchronizing Brain Waves Bring Harmony?

If enhancing memory with light and sound seems futuristic, then welcome to the future. Or so some scientists say. Results from four early-stage clinical trials on mild Alzheimer's disease were presented at the AD/PD 2021 conference, held virtually March 9 to 14. The studies used two closely related approaches to modulate brain waves. Both reportedly synchronized neuronal firing activity in the gamma frequency range, strengthened neuronal connectivity—and perhaps strengthened memory by a bit.

Li-Huei Tsai, Massachusetts Institute of Technology, Boston, and colleagues some years ago turned light and sound into a therapy. They called it GENUS, short for gamma entrainment using sensory stimuli. Gamma brain waves, thought to be a neurophysiological correlate of attention and sensory processing, weaken in people with AD and in mouse models of amyloidosis (Dec 2016 news). Flashing a 40 Hz light at mice for an hour daily for a week enhanced and synchronized the mice's gamma rhythms, which in turn rallied microglia to mop up plaques in their visual cortices. Adding a 40 Hz buzzing sound spread the benefits across the brain, cleared plaques and tangles, and improved the mice's memory (Mar 2019 news; May 2019 news).

Other researchers have partially replicated these results. Sylvain Williams, McGill University, Montreal, and colleagues saw restored hippocampal gamma waves and slightly better spatial memory in amyloidosis mice after 40 Hz flashes using optogenetic stimulation, but plaques did not budge (Etter et al., 2019). Researchers led by Shuzo Sakata, University of Strathclyde, Scotland, U.K., actually saw plaques grow, also after optogenetic stimulation, but they did not assess changes in memory (Wilson et al, 2020).

While a one-time treatment of either light or sound alone boosted people's brain waves, the combination engaged the most brain areas in all participants. The flashing lights and humming sounds triggered no headaches, vision or hearing changes, or seizures, Tsai reported.

Brain Waves May Help Diagnose Dementia

A recently developed two-minute EEG test effectively picks up on early memory impairments associated with Alzheimer's disease.

Flick through the image gallery on your phone or fast forward through a previously watched movie, and notice how the briefly presented images trigger thoughts and memories. This passive visual memory system has been the basis for the design of a brainwave test that may help diagnose Alzheimer's disease.

Alzheimer's & EEG Brain Wave Testing

An EEG is used to measure brain activity while a patient views a series of images of everyday objects on a computer screen. The patient then watches a stream of different images, periodically interspersed with one of the images they saw initially. The images are presented on the screen at a rate of three images per second.

The pictures patients were shown first appear as every fifth image as the test progresses. If the person remembers the image, the EEG readout shows a distinct neural response, and another one five images later, and so on. By looking at the EEG signals for activity at this specific frequency (0.6Hz), we can measure the strength of a person's memory response to previously seen images.

The beauty of this test is that it measures memory passively and objectively. People don't provide any response and don't even need to understand the task. They simply watch a series of pictures while wearing an EEG cap.

Quick & Non-Invasive Testing for Dementia

The EEG procedure takes about 2 minutes. The short time makes the test ideal for use with people with cognitive impairment who can't follow complex task instructions.

By the time a person is diagnosed with Alzheimer's disease – which is the most common form of dementia – they have usually had the disease for years. Diagnosis is not occurring early enough—sometimes a person will have had Alzheimer's for 20 years before it is diagnosed!

Alzheimer's is diagnosed using a combination of reports of memory decline from the patient and their family or carers, often involving tests administered in a clinic. These tests aren't ideal because the patient's anxiety can interfere with the result. They also require the person to be able to speak and write, which makes them ineffective for some people.

Physical Fitness in Institutionalized Older Adults with Dementia: Association with Cognition, Functional Capacity and Quality of Life

Physical fitness, including strength, flexibility, agility/dynamic balance, and aerobic endurance, is important for cognitive function, functional capacity, and quality of life in institutionalized older adults with dementia, suggesting that exercise-based therapeutic strategies should be implemented to improve physical fitness.

Arnaldina Sampaio and colleagues (2020) described physical fitness in institutionalized older adults with dementia. This study examined the association between physical fitness components and cognitive function, functional capacity, and quality of life in institutionalized older adults with dementia. The results showed that all physical fitness components were significantly associated with quality of life. Body mass index was positively associated with caregivers' perception and global quality of life. Aerobic endurance was found to be the key factor in physical fitness associated with cognitive function, functional capacity, and quality of life. The study suggests that strategies to preserve physical fitness components may help maintain cognitive function, functional capacity, and quality of life in this population.

102 Participants were included in the research.

Aspects of the researchers' conclusions appear to back up previous studies in this area: "Agility/dynamic balance is positively associated with global quality of life, according to a multivariate linear regression. This supports previous research linking mobility and quality of life in older adults without dementia," Sampaio suggested.

The researchers suggest that further research is needed to develop efficient exercise programs for institutionalized older adults with dementia. The study aimed to clarify the association between physical fitness components and cognitive function, functional capacity, and quality of life in this population.

OBJECTIVES OF THE TRAINING PROTOCOL

Ensure a Safe and Comfortable Training Environment:

- Establish a secure and inviting training area where participants can practice Table Tennis skills without concerns for their safety.
- Regularly assess the training environment for hazards and maintain a comfortable atmosphere.

Promote Cognitive Functioning through the Table Tennis activities

- Develop Table Tennis training sessions that incorporate cognitive challenges, encouraging strategic thinking, focus, and quick decision-making.
- Create a stimulating and mentally engaging training curriculum that enhances cognitive abilities.

Assist Participants with Activities of Daily Living (ADLs):

- Provide personalized assistance to participants for their daily living activities, such as dressing, grooming, and eating.
- Integrate Table tennis Protocol techniques into daily activities to improve physical coordination and independence.

Deliver Patient Education and Support

- Offer comprehensive education about Table Tennis, including its rules, techniques, and physical benefits.
- Provide continuous coaching and emotional support to participants, helping them overcome challenges and celebrate progress.

Implement Medical Management and Monitoring:

- Collaborate closely with healthcare professionals to assess participants' health status and address any medical needs.
- Monitor and administer necessary medical treatments, medications, or therapies, ensuring participant safety during Table Tennis training.

Provide Caregiver Relief Through TT Protocol Sessions:

- Organize regular Table Tennis Protocol training sessions that offer respite to caregivers, allowing them to take a break while knowing their loved ones are engaged in meaningful activities.
- Establish a flexible training schedule to accommodate caregivers' needs and provide them with the opportunity to recharge.

REQUIREMENTS FOR TABLE TENNIS TRAINING PROTOCOL

1. Qualified Instructors:

- Instructors should be trained in the use of PISTA devices and integrate them effectively into the training sessions.

2. Dedicated Training Space:

- A well-equipped and designated training area with sufficient table tennis tables, nets, and playing surfaces.
- Adequate lighting and ventilation to ensure a comfortable environment for training.

3. Training Curriculum:

- Development of specialized training modules that leverage data for performance analysis and improvement.
- Integration of data to provide real-time feedback on technique and gameplay.

4. Individualized Coaching:

- Utilization of data to tailor coaching and feedback to the specific needs and skill levels of each participant.
- Regular analysis of data to assess individual progress.

5. Regular Practice Sessions:

- Scheduled practice sessions that accommodate participants' availability and skill levels.
- A mix of individual and group training sessions for varied learning experiences.

6. Medical Support and Monitoring:

- Availability of trained medical personnel or access to emergency medical services during training sessions.
- Monitoring of participants' health and physical conditions to prevent overexertion or health issues during training.

7. Record Keeping and Data Analysis:

- Comprehensive data storage and analysis systems to track participant improvements over time.
- Integration of data analysis results into individualized training plans.

8. Caregiver Involvement:

- Involvement of caregivers in PISTA data review sessions, helping them understand and support participants' progress.

9. Evaluation and Feedback:

- Regular evaluation and feedback sessions to gather input from participants and caregivers for continuous program improvement.

PISTA Ping Pong Guidelines

1. Ping Pong without PISTA Device (One Group)

Warm-Up: Start with a gentle warm-up to prepare their bodies. Simple arm stretches and wrist rotations are excellent.

Basic Rules: Depending on the cognitive abilities of the patients, introduce basic Ping Pong rules. Explain that the goal is to hit the ball over the net and into the opponent's side of the table.

Ping Pong Movements:

- **Forehand Stroke:** Teach them how to use a forehand stroke, swinging the paddle forward when hitting the ball. Encourage them to keep their eye on the ball.
- **Backhand Stroke:** Similarly, show them the backhand stroke, where the paddle swings backward when hitting the ball.
- **Serving:** Depending on their abilities, introduce serving techniques. For simplicity, you can have them practice a gentle underhand serve.

Routines: Create routines that include different Ping Pong movements. Here are some examples:

- **Ping Pong Volley:** Patients stand on opposite sides of the table and take turns hitting the ball back and forth, practicing their forehand and backhand strokes.
- **Ping Pong Target Practice:** Set up targets on the table, like plastic cups, and have patients aim to hit the ball into the cups. This adds an element of fun and precision.
- **Ping Pong Relay:** Patients can form teams and take turns hitting the ball over the net, trying to keep a rally going. This encourages teamwork and social interaction.
- **Ping Pong and Memory:** Use colored Ping Pong balls and ask patients to call out the color of the ball they hit or the direction they want to hit (left, right, center) to engage their cognitive skills.
- **Cool Down:** After the Ping Pong session, engage in a cool-down period with gentle stretches to prevent muscle soreness.

Adjust as Needed: Be flexible and adapt the routine to each individual's abilities. Some may need more assistance, while others may be more independent.

2. Ping Pong with PISTA Device (One Group)

Warm-Up: Start with a gentle warm-up, including stretches and wrist rotations, to prepare the participants.

Basic Rules: Explain the basic rules of Ping Pong as mentioned earlier.

Ping Pong Movements: Teach participants the Ping Pong movements, including forehand and backhand strokes.

PISTA Device Integration:

- **Auditory Cues (Beeps):** Program the PISTA device to emit beeping sounds through the headphones. These beeps can serve as cues for when to hit the ball or indicate the start of a rally.
- **Tactile Feedback (Vibration Capsules):** Attach vibration capsules to the participants' wrists or clothing. Configure the PISTA device to trigger vibrations when they successfully make contact with the ball.

Routines with PISTA:

- **Ping Pong Volley with PISTA:** Participants stand on opposite sides of the table, and the PISTA device provides auditory cues (beeps) for when to hit the ball. Simultaneously, vibrations from the capsules reinforce successful hits, providing a multisensory experience.
- **Ping Pong Target Practice with PISTA:** Set up targets on the table, like plastic cups. The PISTA device can beep to guide participants' timing, and vibrations can indicate accurate hits on the target.

Cool Down: After the Ping Pong session, engage in a cool-down period with gentle stretches.

Social Interaction: Encourage conversation and social interaction during breaks.

Scoring with PISTA: The PISTA device can also be programmed to provide different types of beeps or varying vibration intensities to indicate different levels of success. For example, a higher-pitched beep or stronger vibration could signify a particularly good hit.

3. No Table Tennis (one group)

The following guidelines outline a specialized PISTA training program for Alzheimer's patients, excluding the use of table tennis. This program is designed to improve mobility, cognitive functions, and emotional well-being. It focuses on the utilization of the PISTA device, which combines auditory stimulation with vibrations.

Preparation:

- **Assessment and Personalization:** Start with a thorough assessment of each patient's condition, including mobility, cognitive abilities, and emotional state. Tailor the PISTA program to their specific needs and goals.
- **Consultation:** Consult with licensed PISTA coach experienced in PISTA training to ensure that the program aligns with the patient's overall treatment plan.

Training Routine:

- **Regular Schedule:** Establish a consistent training routine, incorporating PISTA sessions into the patient's daily activities.
- **Device Setup:** Ensure that the PISTA device is set up correctly, with both earphones and vibration capsules properly connected. Adjust settings according to the patient's preferences.
- **Daily Session:** Incorporate PISTA sessions as part of the daily routine. These sessions should be scheduled for a duration that suits the patient's stamina and preference.

Physical Mobility:

- **Gait and Movement:** During PISTA sessions, integrate walking or seated mobility exercises. The combination of auditory stimulation and vibrations can enhance coordination and mobility.
- **Stretching:** Include gentle stretching exercises to improve flexibility and reduce muscle stiffness. These exercises can help patients feel more comfortable in their daily activities.

Cognitive Stimulation:

- **Memory and Cognitive Exercises:** Engage in PISTA programs designed to enhance cognitive functions. These exercises can include memory games, problem-solving tasks, and activities to boost attention and focus.
- **Mindfulness and Relaxation:** Integrate meditation and mindfulness programs to encourage mental relaxation and emotional well-being.

Emotional Well-being:

- **Stress Reduction:** Utilize PISTA's relaxation programs to help patients manage stress and anxiety. Regular practice can contribute to improved emotional well-being.
- **Motivation and Positivity:** Leverage PISTA to promote a positive outlook. Encourage patients to set personal goals and use the device's motivational programs to stay inspired.
- **PISTA Modalities:** Use of different modalities such as P15, P5 and Entry Points

Monitoring and Progress:

- **Data Collection:** Track the patient's progress by recording changes in mobility, cognitive abilities, and emotional state.
- **Regular Assessments:** Conduct regular assessments with the help of a PISTA coach to evaluate the effectiveness of the PISTA program.
- **Feedback and Adjustments:** Based on assessments and patient feedback, make necessary adjustments to the program, such as changing exercises, durations, or settings.

Caregiver Involvement:

- **Training and Supervision:** Caregivers should be trained on PISTA device usage to assist patients. They can help set up the device and guide patients through exercises.
- **Motivational Support:** Caregivers play a crucial role in providing motivation and emotional support, encouraging patients to adhere to their training regimen.

4. Ping Pong with PISTA Device – Home use (One Group)

The following guidelines are designed for Alzheimer's patients who want to combine table tennis and PISTA device training for use at home. This program aims to enhance mobility, cognitive function, and emotional well-being. Table tennis and the PISTA device will be integrated to provide a comprehensive training experience.

Preparation:

- **Assessment and Personalization:** Start with a thorough assessment of each patient's condition, including mobility, cognitive abilities, and emotional state. Tailor the PISTA program to their specific needs and goals.
- **Consultation:** Consult with licensed PISTA coach experienced in PISTA training to ensure that the program aligns with the patient's overall treatment plan.

Training Routine:

- **Regular Schedule:** Establish a consistent training routine, incorporating table tennis and PISTA sessions into the patient's daily activities at home.
- **Device Setup:** Ensure that the PISTA device is set up correctly, with both earphones and vibration capsules properly connected. Adjust settings according to the patient's preferences.
- **Table Tennis Equipment:** Provide access to a table tennis setup, including a table, paddles, and balls. Ensure the area is safe for playing.

Physical Mobility and Table Tennis:

- **Table Tennis Sessions:** Incorporate regular table tennis sessions as part of the training program. Table tennis is an excellent way to improve hand-eye coordination, reflexes, and overall mobility.
- **Exercise Variety:** Include table tennis drills and practice exercises that mimic the demands of the game. Patients can work on their forehand and backhand strokes, footwork, and rallying skills.
- **Combining PISTA:** For table tennis sessions, use the PISTA device's sound and vibration to improve focus and coordination during the game. Adjust the PISTA settings for optimal enhancement.

Monitoring and Progress:

- **Data Collection:** Track the patient's progress by recording changes in mobility, cognitive abilities, and emotional state.
- **Regular Assessments:** Conduct regular assessments with the help of a PISTA Coach to evaluate the effectiveness of the combined program.
- **Feedback and Adjustments:** Based on assessments and patient feedback, make necessary adjustments to the program, such as changing exercises, durations, or settings.

Caregiver Involvement:

- **Training and Supervision:** Caregivers should be trained on both the PISTA device and basic table tennis rules and techniques. They can assist patients in setting up the PISTA device and ensure that table tennis sessions are conducted safely.
- **Motivational Support:** Caregivers play a crucial role in providing motivation and emotional support, encouraging patients to adhere to their training regimen.

TABLE TENNIS TRAINING SESSIONS ROUTINE AND INSTRUCTIONS

Notes:

- Throughout the training session, participants should try to replicate the movements from the videos, with eyes closed if possible.
- Ensure that all key points from the text are incorporated into the training.
- Focus on maintaining good posture and balance during the walking forward test.
- Regularly review and compare the filmed sessions to track progress.
- Adapt the training plan as needed based on test results and performance feedback.

Using EEG Device before training

Before starting table tennis training, the patient will use an EEG device to measure their brain activity to understand its condition before the session. They will repeat this during breaks after the initial training phase and at the end of the entire session. We will record these brainwave measurements throughout. Patients don't need to respond or understand the task, and the EEG procedure only takes about 2 minutes.

1. Warm-Up (5 minutes):

Ensuring that the trainees have an effective warm-up before each table tennis training session is essential for their safety, performance, and overall development. Here are some instructions and guidelines to help you lead your trainees through a productive 5-minute warm-up routine:

a. Importance of Warm-Up:

Emphasize the significance of the warm-up to prevent injuries and enhance performance.

b. Clear Communication:

Clearly explain each warm-up exercise and demonstrate proper technique.

c. Progressive Intensity:

Begin with low-intensity activities, gradually increasing to moderate intensity in the final minute.

d. Mental Preparation:

Encourage trainees to visualize their game and set goals during the warm-up.

2. First Session (1 hour):

a. Filming Players Walking Forward and Analyzing Posture

- Begin by explaining the importance of posture in table tennis.
- Film each trainee as they walk forward. Focus on their posture, balance, and body alignment.
- Provide immediate feedback on posture, highlighting areas for improvement.

b. Footwork Training

- Gather trainees to watch instructional videos featuring Coach Xu Gang and Pista Trainer.
- Pay attention to the footwork techniques demonstrated in the videos.
- Encourage trainees to take mental notes and ask questions for clarification.

c. Basic Training Exercises

- Transition to practical training exercises based on the video demonstrations.
- Start with basic footwork drills, ensuring that trainees implement the techniques they observed.
- Emphasize proper form and coordination.

d. Filming the Training Session

- Pause briefly to film a segment of the training session in progress.
- Explain that this footage will be used for later analysis to track progress.

e. Repeat Walking Forward Test and Evaluate Balance

- Have trainees repeat the walking forward test.
- Assess whether there have been any improvements in posture and balance.

f. Conduct the Sit-Down-Stand-Up Test

- Initiate the sit-down-stand-up test to evaluate lower body strength, balance, and agility.
- Observe trainees' performance and note any challenges or improvements.

g. Recap and Closing

- Summarize the key takeaways from the session, focusing on posture, footwork, and balance.
- Encourage trainees to practice what they've learned between sessions.
- Highlight the importance of consistency in training.

3. Break (15 minutes):

Rest and rehydrate.

4. Second Session (30 minutes):

a. Continue Working on Footwork and Posture

- Begin with a brief recap of the importance of footwork and posture in table tennis.
- Conduct footwork drills to reinforce the techniques learned in the first session.
- Emphasize proper posture, balance, and quick movement.

b. Play with the Coach in Returning the Ball

- Engage in practice rallies with the trainees, acting as their opponent.
- Focus on returning balls with different spins, speeds, and angles to challenge the trainees' footwork and stroke techniques.
- Provide immediate feedback on their shot placement and footwork.

c. Filming the End of the Session

- Pause briefly to film a segment of the training session's closing phase.
- Explain that this footage will be used for further analysis and feedback.

d. Recap and Closing

- Summarize the key points covered in the session, highlighting improvements in footwork and posture.
- Encourage trainees to continue practicing and refining their skills.
- Discuss any specific areas for further improvement or focus in the next session.

5. Post-Training Session Analysis

a. Review the Recorded Training Session:

- Start by reviewing the video footage of the entire training session. Pay close attention to each trainee's performance, posture, footwork, and overall technique.

b. Assess Progress and Areas for Improvement:

- Evaluate each trainee's progress compared to previous sessions and goals set.

- Identify specific areas where trainees have shown improvement and areas that require further development.
- Take note of any consistent patterns or challenges among the trainees.

c. Individual Assessment:

- Provide individualized feedback for each trainee based on their performance during the session.
- Highlight strengths and areas requiring improvement.
- Address any concerns or questions raised by the trainees during the session.

d. Set Goals for Future Training Sessions:

- Based on your analysis, establish specific training goals for the next session.
- Prioritize areas that require additional attention and set objectives to work on those aspects.

e. Adjust Training Plans:

- Modify training plans, drills, or exercises to better address the identified areas of improvement.
- Ensure that future sessions are tailored to meet the evolving needs of the trainees.

f. Communicate with Trainees:

- During the next training session, share the results of the analysis with the trainees.
- Encourage open discussion and feedback regarding their performance and training experiences.

g. Document the Analysis:

- Maintain a record of your analysis, including individual progress notes and collective observations.
- Use this documentation as a reference point for tracking trainee development over time.

6. Test

The objective of this test is to measure the progress and improvements made by our trainees over the course of the training program. We will focus on key areas, including posture, footwork, ball control, and overall performance.

Test Components:

- **Posture Assessment:** Evaluate each trainee's posture, body positioning, and balance during play. Note any improvements compared to their initial assessments.
- **Footwork Evaluation:** Assess the trainees' footwork techniques, agility, and movement around the table. Consider their ability to position themselves for optimal shots.
- **Ball Control Test:** Observe the trainees' control over the ball, including their ability to return different types of shots (e.g., spins, speed variations) with accuracy.
- **Match Play:** Organize practice matches among trainees to assess their overall game strategy, shot selection, and decision-making under real game conditions.

Scoring and Feedback:

- **Use a Scoring System:** Assign scores or ratings for each aspect assessed, such as posture, footwork, ball control, and match performance.
- **Provide Feedback:** After the test, offer individualized feedback to each trainee based on their performance. Highlight areas of improvement and acknowledge their strengths.

Documentation:

- **Record Results:** Keep detailed records of each trainee's performance during the test, including scores, observations, and feedback.
- **Compare to Previous Assessments:** Compare the results of this test to the initial assessments conducted at the beginning of the training program. This will help measure progress effectively.

TABLE TENNIS TRAINING VIDEOS

- Dr Andre Stang Training Videos
- PISTA Ping Pong Footwork movements
- PISTA Ping Pong for Athletes
- PISTA Ping Pong for Alzheimer's Patients.

ALL TRAINING VIDEOS FOR TABLE TENNIS WILL BE SENT THROUGH WHATSAPP



PISTA MODALITIES

The use of various modalities within the PISTA device, including P15, P25, P5, and the strategic utilization of entry points, plays a pivotal role in the management of emotions, including fears and anxiety, for patients with Alzheimer's. P15 modalities offer a tailored approach, addressing specific emotional objectives and enabling patients to experience transformative personal growth. These customizable exercises empower patients to navigate their emotions, such as sadness and anxiety, with a newfound understanding, ultimately fostering emotional healing and a renewed sense of well-being.

The integration of PISTA with these modalities enhances their effectiveness. PISTA's auditory beat stimulation and vibrations create a sensory-rich environment that engages individuals more deeply with these exercises. By synchronizing the auditory and vibrations with P15 modalities, patients experience sensory feedback that enhances emotional exploration and regulation. PISTA acts as an enabler, making these exercises more immersive and impactful.

The P5 + 7 motivation technique synergizes seamlessly with PISTA. This technique promotes a sense of motivation, inspiring individuals to overcome their emotional challenges, including fears and anxiety. PISTA complements this by introducing auditory and vibration stimuli synchronized with the P5 + 7 technique. This multi-sensory approach provides a heightened level of motivation, encouraging individuals to move from planning to action with determination.

Entry points, integrated with PISTA, offer structured anchors for emotional management. With the addition of synchronized auditory and vibration from PISTA, these entry points become even more effective. Patients can address their emotional challenges with more focused support, and the sensory integration enhances the overall experience of working through fears and anxieties.

In summary, the PISTA device enhances the emotional management modalities by providing synchronized sensory experiences. This integration ensures that the emotional exercises are more engaging, motivating, and effective for individuals with Alzheimer's. By leveraging PISTA alongside these modalities, patients can address their emotional well-being more comprehensively, leading to a renewed sense of hope and an improved quality of life.

1. P15

15 Minutes Use of the PISTA Machine

Each P15 application is designed to achieve specific and targeted results. By tailoring the P15 modalities to address particular objectives, the patient can effectively work towards their desired outcomes. Through the implementation of these customized P15 exercises, the patient gains valuable insights and experiences transformative changes, enabling them to make meaningful progress in their personal and emotional growth. The flexibility and adaptability of the P15 modalities ensure that the patient's unique needs and goals are addressed, leading to a more fulfilling and successful journey of self-improvement.

By implementing the P15 techniques, the patient's brain is expected to undergo training to process sadness, gaining a deeper understanding of its purpose and learning to work through it without becoming emotionally trapped. This process aims to achieve a balance between feeling sad and enjoying daily living, enabling the patient to effectively deal with grief and establish a new, positive relationship with herself and her future self. Ultimately, the program endeavors to bring relief from the burden of sadness and pain, fostering emotional healing and a renewed sense of well-being for the patient.

Step 1.

While having the machine on, the coach will ask the client to close their eyes during the entire 15 minutes of wearing the machine.

Step 2.

The coach will suggest a topic to the client. Tell the client to allow their brain and mind to just focus on the topic.

The topic chosen will be called as the **THEME**.

Sample Topics: The topic that causes them anxiety, stress or a topic that they want a solution or a topic that they want some more insights on or a habit. E.g., Planning your day

Step 3.

Tell the client that don't push himself/herself to think of answers.

If their mind drifts away from the given topic, it's OKAY.

Just allow these 15 minutes to focus on their thoughts.

Step 4.

After 15 minutes, tell them to write down what comes out and what other thoughts after using P15.

Step 5.

Tell them that they can record what topic has done, what are the results, how it helps them and how their emotion changes, and how are they feeling after.

2. P25

- 25 Minutes Use of The PISTA Machine but Should Be Focused Only on A Certain Task
- You have to be away from your mobile or any kind of disturbance you just need to focus on the topic or task or some assignment you need to work on.
- Try to estimate the time needed for this task, choose something that requires 2 hours or three and you can see that you can reduce a lot of time.
- P25 only need 25 minutes of using the PISTA device while working on your certain task but during that period of time you need to be away from any of the other media or computer and you will see how you will benefit from this 25-minute of focus and concentration and get your work done.

Integrating the P25 modality into table tennis training for patients with Alzheimer's represents a targeted strategy to enhance their cognitive functions, concentration, and skill development. To tailor this integration, it is crucial to set clear training objectives that consider the unique challenges associated with Alzheimer's, such as memory retention and cognitive decline. Dedicate a specific 25-minute segment during the training session for the use of the PISTA device with the P25 modality, ensuring a conducive environment with minimal distractions. Select a table tennis training task that aligns with the patient's objectives and effectively addresses the cognitive challenges posed by Alzheimer's.

The PISTA device, alongside the P25 modality, can enhance the training experience by promoting focused engagement in the chosen task. It achieves this by providing synchronized auditory and tactile stimuli, which help patients stay immersed in the training and reduce external distractions. Following the 25-minute session, assess the patient's progress and gather feedback on their experience.

As patients adapt to the exercises, consider gradually increasing the complexity and cognitive demands of the tasks, using the PISTA Device sounds and beeps to challenge different aspects of their Alzheimer's-related cognitive decline. This tailored approach empowers individuals with Alzheimer's to make significant improvements in their table tennis skills while effectively managing the cognitive demands associated with their condition. It fosters a supportive and structured training environment that aligns with their unique needs and goals.

3. P5+7

In the P5 + 7 motivation technique, the patient takes five minutes to carefully contemplate and plan their upcoming activity while utilizing the vibration and sound features at the fastest speed on the Pista device. Once the five minutes have elapsed, the patient is prompted to immediately put their planned action into motion without delay. This modality aims to enhance motivation and prompt swift action, enabling the patient to efficiently transition from planning to execution and achieve their desired goals with focused determination.

- Patients can choose specific table tennis videos to watch, focusing on aspects of the game they aim to improve, such as stroke techniques, footwork, or strategy.
- Begin with the five-minute planning phase, during which patients use the PISTA device's vibration and sound features to concentrate and set their goals. This phase is ideal for taking notes, visualizing strategies, and mentally preparing for what they will learn from the video.
- Patients watch the selected table tennis video, paying close attention to the techniques and insights presented.
- Immediately after watching the video, patients transition to the action phase. They put the insights gained from the video into practice on the table tennis court.
- Patients repeat this process regularly, using different videos to target various aspects of their game. By consistently applying the P5 + 7 modality, they can efficiently learn and improve their table tennis skills while staying motivated and determined.

4. Entry Point

- 11 SETS
- Each set is equivalent to one and a half minutes

Step 1.

While having the machine on, the coach will ask the client to close their eyes during each set. You can tell them to time themselves using a mobile timer.

Step 2.

The coach will give an ENTRY POINT topic to the client. Tell the client to allow their brain and mind to just focus on the topic on each set.

WHAT IS ENTRY POINT?

An entry point is a thought or an image to which the client focuses on throughout the session. The entry point can also be a scenario or an experience in the past, which an individual continues to feel strongly about. The coach helps the client in identifying this image and later in applying the stimulus that will stimulate the brain. The entry point may change as the sessions continue, depending on the progress of the user.

Step 3.

After each set of one and a half minute, tell the client to write down what are their thoughts.

Aside from recording the thoughts, tell them to also record their feelings, emotions or physical emotions that might arise during each set. E.g.: headache, sobbing, weeping laughter, difficulty in breathing etc.

Step 4.

After finishing 11 sets. Tell the client that they need to take a short quiet break of five to ten minutes and let their mind sink in.

Step 5.

Then, the client has to write down the summary of the new insights n feeling they have for the session.

Research Evidence for Management Plan

Exercise Therapy

BALANCE TRAINING

Falls are a major health concern in elderly populations. There is a strong correlation between fall occurrence and balance, which causes notable concern for people with AD as balance issues can appear early in the diagnosis. Significant deficits in both static and dynamic balance can be noted, with specific difficulties in balance during turning, gait and dual task activities. Both individual and group exercise sessions were found to provide benefit. Specific exercises aimed at people with mild to moderate AD include tandem walks, with and without eye tracking, high side steps and weight shifts. Video game interventions also exist that are likely to play a larger role for future generations that develop AD. Specifically, using a Wii-fit intervention was shown to decrease fear of falls and increase Berg Balance Scale scores compared to a walking program. An important criteria of balance programs and exercise programs in general is to be adequately challenging to show noticeable improvements. Therefore, clinicians should continue this protocol even in patients with AD.

STRENGTH TRAINING

Strength Training promotes significant improvements in patients with AD. Studies have shown the improvements were with respect to lower limb strength and balance. Resistance training with ADs not only increases muscle strength by 3-17%, but it also improves ambulatory abilities by increasing gait speed and enhances the individual's ability to perform ADLs. By getting individuals with ADs to perform resistance exercises, we decrease the amount of muscle loss which decreases inactivity and the loss of functional ability. One study showed that with a 60-minute resistance training session 2x a week, balance and muscular endurance improved significantly. Cardiorespiratory improvements were also noted in individuals who performed resistance training which would further promote an active lifestyle in individuals with AD.

Exercises that focus on functional tasks demonstrate a medium effect size on slowing the decline of cognitive ability, improving orientation, memory function and performance of functional task in AD patients. These strategies focus on learning compensatory techniques that can be applied to activities of daily living.

AEROBIC TRAINING

Aerobic Training has been shown to improve an AD's patient's ability to independently perform activities of daily living (ADLs) by improving their overall functional capacity. Studies have also shown that improvements in cardiorespiratory fitness were related to better memory performances. Although currently still unsure of the optimal modality and dosage of aerobic exercise, there is evidence that shows cognitive and physical improvements can occur to help delay the progression of ADs while allowing the patient to live an active life.

A Reducing Disability in Alzheimer Disease (RDAD) program was designed to reduce impairments in independence and behavior in AD populations. The RDAD has an exercise component that consists of 30min/day of aerobic activities, strength training, balance, and flexibility training. A study that compared the RDAD to normal routine care showed improved physical functioning after 3 months in RDAD patients that exercised at least 60min/wk. This effect lasted 2 years after the study. These patients also were less institutionalized.

Balance Exercises

1. Tandem Stance (can progress to Tandem Walk)
2. Weight Shift Side to Side

Frequency: Daily


Intensity: As tolerated

Duration: 1 set of 3 repetitions, 10 seconds hold each side. Progress as needed



Tandem Stance

Stand one foot directly in front of the other and hold for 10 seconds each side or as long as you can. Repeat with other foot in front.



Weight Shift Side to Side

Gently sway (lean) to the right and left so that your weight shifts from your right foot to your left foot. Do not lift your toes or heels. Be sure that your shoulders and hips move together. Do not bend at your hips. Slowly increase how far you can sway right and left without taking a step.

Step Ups

Frequency: Daily

Intensity: As tolerated

Duration: 1 set of 10 repetitions. Progress as needed



Sit to Stands

Frequency: Daily

Intensity: As tolerated

Duration: 1 set of 10 repetitions. Progress as needed



Glute Bridges

Frequency: Daily

Intensity: As tolerated

Duration: 1 set of 10 repetitions, hold for 3 seconds at the top. Progress as needed



Walking

Frequency: Daily

Intensity: Light

Duration: Start at 20 min/day. Progress 5 mins per week or as tolerated

Pushing Exercises

- Wall Push Ups
- Chair Dips

Frequency: 3x/week

Intensity: Moderate

Duration: 1 set of 10 repetitions. Progress as needed



Wall Push Ups

Stand in front of a wall, about two feet away, and place your hands against it at shoulder height. Keeping your body straight, bend at the elbows and lean towards the wall, and then straighten your arms and push back into an upright position.



Chair Dips

Shoulder Elevation

- Lifting objects (ex. canned foods, cups, dishes) from table to overhead
- Resistance band Shoulder Press

Frequency: 3x/week

Intensity: Moderate

Duration: 1 set of 10 repetitions. Progress as needed



Lifting objects



Resistance band Shoulder Press

Neck Mobility Exercise for Patients with Dementia

A. Left and Right

1. Begin by gently guiding the patient's head to the left. Move slowly and patiently.
2. Allow the patient a moment to feel the gentle stretch in their neck as they look to the left.
3. Now, gently guide their head to the right, assisting them in looking over their right shoulder.
4. Give the patient time to experience the gentle stretch as they look to the right.
5. Continue this gentle left-to-right movement alternately for a few repetitions. Adjust the pace to match the patient's comfort level.



Neck Mobility Exercise for Patients with Dementia

B. Looking Down

1. Gently guide the patient's chin downward toward their chest.
2. Encourage the patient to look downward, focusing on their chest or lap. Emphasize that they should keep their chin tucked gently, without tilting their head too far back.
3. After a moment, slowly guide their head back to a neutral, straight-ahead position.
4. Continue this gentle chin tuck movement, alternating between looking downward and returning to a neutral position. Perform a few repetitions at a pace comfortable for the patient.



C. Neck Tilt Exercise

1. Begin by gently guiding the patient's head to the left, tilting their ear toward their left shoulder
2. Now, gently guide their head to the right, tilting their ear toward their right shoulder.
3. Continue this gentle left-to-right neck tilting movement alternately. Ensure a slow and unhurried pace, and adjust according to the patient's comfort level.



Berg Balance Scale (BBS)

The Berg Balance Scale (BBS) is used to objectively determine a patient's ability (or inability) to safely balance during a series of predetermined tasks. It is a 14 item list with each item consisting of a five-point ordinal scale ranging from 0 to 4, with 0 indicating the lowest level of function and 4 the highest level of function and takes approximately 20 minutes to complete. It does not include the assessment of gait.

Equipment required

- A ruler
- 2 standard chairs (one with arm rests, one without)
- A footstool or step
- 15 ft walkway
- Stopwatch or wristwatch
- PISTA Device

Name: _____ Date: _____

ITEM DESCRIPTION SCORE (0-4)

1. Sitting to standing _____
2. Standing unsupported _____
1. Sitting unsupported _____
2. Standing to sitting _____
3. Transfers _____
4. Standing with eyes closed _____
5. Standing with feet together _____
6. Reaching forward with outstretched arm _____
7. Retrieving object from floor _____
8. Turning to look behind _____
9. Turning 360 degrees _____
10. Placing alternate foot on stool _____
11. Standing with one foot in front _____
12. Standing on one foot _____

Total _____

Berg Balance Scale (BBS)

General instructions for completing the scale

Please document each task and/or give instructions as written. When scoring, please record the lowest response category that applies for each item.

In most items, the subject is asked to maintain a given position for a specific time. Progressively more points are deducted if:

- the time or distance requirements are not met.
- the subject's performance warrants supervision.
- the subject touches an external support or receives assistance from the examiner.

The subject should understand that they must maintain their balance while attempting the tasks. The choices of which leg to stand on or how far to reach are left to the subject. Poor judgment will adversely influence the performance and the scoring.

Equipment required for testing is a stopwatch or watch with a second hand, and a ruler or other indicator of 2, 5, and 10 inches. Chairs used during testing should be a reasonable height. Either a step or a stool of average step height may be used for item # 12.

Interpretation

Cut-off scores for the elderly were reported by Berg et al 1992 as follows:

- A score of 56 indicates functional balance.
- A score of < 45 indicates individuals may be at greater risk of falling.

It has been reported more recently that in the elderly population a change of 4 points is needed to be 95% confident that true change has occurred if a patient scores within 45–56 initially, 5 points if they score within 35–44, 7 points if they score within 25–34 and, finally, 5 points if their initial score is within 0–24 on the Berg Balance Scale.

FIVE TIMES SIT TO STAND TEST (FTSTS)

Description: Assesses functional lower extremity strength, transitional movements, balance, and fall risk.

Equipment: Stopwatch; standard height chair with straight back (16 inches high);

Therapist Instructions: Have the patient sit with their back against the back of the chair. Count each stand aloud so that the patient remains oriented. Stop the test when the patient achieves the standing position on the 5th repetition.

Patient Instructions: "Please stand up straight as quickly as you can 5 times, without stopping in between. Keep your arms folded across your chest. I'll be timing you with a stopwatch. Ready, begin."

Auditory Cue (Beeps): As the patient starts performing the sit-to-stand repetitions, the PISTA Device should emit a series of beeps through the headphones. These beeps can serve as an auditory cue to help patients maintain their pace and rhythm during the test. The beeps should be timed to match the desired speed of the test.

Vibrational Feedback (Vibration Capsules): In addition to the auditory cues, the PISTA Device should also incorporate vibration capsules strategically placed on the patient's body, such as on their thighs or back. These capsules should vibrate at key points during the test, providing tactile feedback to assist the patient in maintaining balance and control during the sit to-stand movements.

Scoring the test

- Record the time from command "Go" until the patient's buttocks touch the seat following the 5th stand.
- Perform 1-3 timed trials.
- Final score is the average of the timed trials.

Interpretation:

- **Lower times = better scores**
- **MDC:** 3.6 - 4.2 secs
- **MCID:** 2.3 secs

Age-Matched Norms

Age Bracket	Time (sec)
60-69 yo	11.4
70-79 yo	12.6
80-89 yo	14.8

TIMED UP & GO (TUG)

Purpose: To assess mobility

Equipment: A stopwatch

Directions: Patients wear their regular footwear and can use a walking aid, if needed. Begin by having the patient sit back in a standard arm chair and identify a line 3 meters, or 10 feet away, on the floor.

1. Instruct the patient:

When I say "Go," I want you to:

- Stand up from the chair.
- Walk to the line on the floor at your normal pace.
- Turn.
- Walk back to the chair at your normal pace.
- Sit down again.

2. On the word "Go," begin timing.

3. Stop timing after patient sits back down.

4. Record time.

Time in Seconds: An older adult who takes ≥ 12 seconds to complete the TUG is at risk for falling.

Observe the patient's postural stability, gait, stride length, and sway.

Check all that apply:

- Slow tentative pace
- Loss of balance
- Short strides
- Little or no arm swing
- Steadying self on walls
- Shuffling
- En bloc turning
- Not using assistive device properly

These changes may signify neurological problems that require further evaluation.

Caregiver Assessment

Caregiver Assessment Protocol

Objective: To conduct an assessment of patients after they have used PISTA and provide v
Procedure:

A. Assessment Preparations:

- **Schedule Assessment:** Coordinate with the PISTA coach to establish a schedule for assessments that aligns with the patient's training sessions.
- **Collect Assessment Tools:** Gather relevant assessment tools, which may include self-report forms, questionnaires, or any specific guidelines provided by the PISTA coach.
- **Patient Briefing:** Before conducting the assessment, communicate with the patient. Explain the purpose of the assessment and assure them that their feedback is valuable for enhancing their experience with PISTA.

B. Conducting the Assessment:

Physical Evaluation:

- Observe the patient's physical state, including their balance, mobility, and coordination.
- Note any improvements or challenges they may face after using PISTA.

Cognitive Assessment:

- Assess the patient's cognitive abilities, such as memory, attention, and problem-solving skills.
- Determine if there have been any noticeable cognitive enhancements.

Dual Tasking Evaluation:

- Engage the patient in dual tasking scenarios to observe their multitasking capabilities and identify any improvements or areas that need further attention.

C. Recording and Reporting:

- **Collect Data:** Compile the patient's self-reports, physical observations, cognitive assessments, and dual tasking results.

- **Analysis:** Analyze the gathered data to identify trends, improvements, or challenges the patient may be experiencing after using PISTA.
- **Consult with the Patient:** Discuss the assessment findings with the patient, clarifying any areas of concern or potential improvements.
- **Complete Assessment Forms:** Fill out the caregiver assessment forms or any reporting templates provided by the PISTA coach.

D. Reporting to PISTA Coach:

- Share the completed assessment forms or reports with the PISTA coach through the designated communication channels.
- Schedule a meeting or discussion to review the assessment findings and provide context.

Post-Assessment Follow-up:

Collaborate with the PISTA coach to develop an action plan based on the assessment results. This plan may involve adjusting PISTA sessions or introducing new techniques.

Continue to engage with the patient, monitor their progress, and provide any additional information or observations to the PISTA coach for ongoing optimization.

Dual Tasking Assessment:

- **Explain Dual Tasking:** Clarify the concept of dual tasking, which involves performing two tasks simultaneously, such as walking and counting, to evaluate cognitive and motor skills.
- **Choose Appropriate Dual Tasks:** Select dual tasks that are relevant to patients' daily lives, such as walking while talking or carrying objects.
- **Demonstrate Dual Tasking:** Demonstrate how to perform the chosen dual tasks. Ensure that patients understand the tasks and their objectives.
- **Start with Simple Tasks:** Begin with simple dual tasks and gradually progress to more complex ones as patients become more comfortable with the process.
- **Safety Precautions:** Emphasize safety during dual tasking. Patients should prioritize their physical safety, and you should monitor them to prevent any accidents.
- **Recording Observations:** Keep detailed notes on patients' performance during dual tasking. Note any challenges, such as balance issues or difficulties in managing both tasks simultaneously.

- **Feedback and Improvement:** After completing dual tasks, provide feedback to patients. Discuss any areas where improvement may be needed and suggest strategies to enhance their dual tasking abilities

PISTA SOUND AND VIBRATION INTEGRATION WITH THE TRAINING PROTOCOLS

The integration of PISTA Sound and Vibration technology, driven by the potential of binaural beats, represents a groundbreaking advancement in designing exercise protocols specifically tailored to meet the unique needs of patients grappling with Alzheimer's and dementia. PISTA technology introduces a distinctive sensory experience by seamlessly intertwining the power of binaural beats with subtle vibrations, ushering in a transformative approach to exercise regimens for these individuals.

Binaural beats, at the heart of this innovative technology, offer a remarkable auditory stimulus that goes beyond mere sound. They possess the ability to heighten patient awareness of their physical movements while playing a pivotal role in supporting motor control. For individuals contending with Alzheimer's and dementia, these beats serve as vital cues, aiding in the initiation and maintenance of specific motions. This is particularly significant as cognitive challenges can often result in difficulties with initiating actions or responding to cues.

Additionally, the profound impact of PISTA technology extends beyond the physical realm. It wields the capacity to provide rhythmic cues that go far in enhancing cognitive functions. This support spans various aspects of cognition, including memory, attention, and overall cognitive well-being. By entraining the brain through binaural beats, patients experience a unique form of cognitive stimulation that fosters mental clarity and engagement.

What sets PISTA technology apart is its ability to offer immediate feedback on the quality of movements, thereby motivating patients to make real-time adjustments. This feedback loop between the technology and the individual is pivotal in maintaining the patient's engagement and fostering a sense of empowerment.

Furthermore, the interactive and engaging nature of PISTA technology, when combined with binaural beats, significantly enhances patient adherence to exercise routines. This is achieved through personalized therapy plans that cater to individual needs and preferences. The data collected during these sessions serves a dual purpose: it aids in the monitoring and assessment of progress while contributing to the refinement of the therapy plan.

Crucially, PISTA technology's seamless integration into a user-friendly device facilitates at-home therapy and remote supervision, ensuring a consistent continuum of care for patients grappling with Alzheimer's and dementia. This democratizes access to rehabilitation and support, irrespective of geographical constraints.

In conclusion, PISTA Sound and Vibration technology, fortified by the incorporation of binaural beats, stands as a beacon of hope for individuals dealing with Alzheimer's and dementia. Beyond its physical benefits, it offers a lifeline of cognitive and emotional support. By actively engaging in this rehabilitation journey, patients experience tangible improvements in cognitive functions and overall quality of life, reaffirming the potential for meaningful progress and enhanced well-being.

BINAURAL SOUND

The integration of binaural sound through PISTA technology holds great promise for enhancing the well-being of individuals with Alzheimer's disease. Binaural sound engages and mentally stimulates patients during exercises, reducing monotony and aiding in stress reduction. It also supports cognitive functions, enhancing focus and concentration, which are often challenging for Alzheimer's patients. Additionally, it uplifts mood and motivation for physical activities. This holistic approach fosters a stronger mind-body connection, improving coordination and body awareness. In summary, binaural sound in training protocols complements Alzheimer's therapy by addressing both mental and physical aspects, ultimately enhancing the overall quality of life.

VIBRATIONS

The integration of binaural beats through PISTA technology offers a holistic approach to enhancing the well-being of individuals with Alzheimer's. It engages auditory and tactile senses, fostering heightened awareness of the body and its movements. This approach aids in initiating and maintaining specific motions, particularly beneficial for individuals with difficulties in starting movements due to Alzheimer's. Additionally, it contributes to balance and coordination training and enhances focus, motivation, and stress reduction. The flexibility of this approach allows for personalized therapy plans, and data collected helps track progress for optimization. In summary, the integration of binaural beats through PISTA technology enhances the physical and cognitive aspects of Alzheimer's therapy, improving overall quality of life.

PISTA DEVICE INTEGRATION FOR COGNITIVE AND PHYSICAL WELL-BEING

Integrating the PISTA device into the daily routines and training exercises of Alzheimer's patients, with the invaluable assistance of their caregivers, offers a comprehensive approach to enhancing well-being. Caregivers can play a pivotal role in facilitating this integration, ensuring that patients receive optimal benefits from the technology.

Each day can begin with caregivers assisting patients in donning the PISTA device, including earphones, as part of their morning routine. This collaborative effort helps patients experience synchronized auditory stimulation, setting a positive tone for the day ahead. During exercise sessions, caregivers can provide guidance, helping patients adjust the device settings, select appropriate programs, and maintain motivation throughout the workout.

In cognitive training exercises, caregivers can actively participate alongside patients, both utilizing the PISTA device to enhance cognitive engagement. This shared experience not only bolsters cognitive functions but also strengthens the emotional connection between patients and caregivers.

For exercises aimed at improving balance and mobility, caregivers can offer physical support while patients wear the PISTA device, ensuring safety and confidence during these activities. Meditation and relaxation sessions can become shared moments of calm, with both caregivers and patients benefiting from the deepened meditative state facilitated by PISTA.

Caregivers can also take an active role in coordinating and personalizing the PISTA device's programs, collaborating with healthcare providers to tailor settings and programs to the patient's evolving needs and goals. They can assist in monitoring and tracking the patient's progress, recording data, and sharing it with healthcare providers for ongoing assessment and adjustments.

Perhaps most importantly, caregivers provide consistent motivation, encouragement, and emotional support throughout the patient's exercise and training routines with the PISTA device. This steadfast support significantly enhances patient motivation and adherence to the regimen.

By integrating the PISTA device into daily life and training exercises, with caregivers as partners in this journey, patients with Alzheimer's can unlock the full potential of this innovative technology. Together, they can enhance cognitive function, emotional well-being, and overall

quality of life, while simultaneously strengthening their bond through shared experiences and mutual support.

HOW THE PISTA METHOD WORKS

The PISTA methodology employs a range of techniques, harnessing the power of brainwave feedback, brain stimulation, and self-regulation procedures to optimize brain function and correct inefficient mental activities, with a specific focus on addressing the unique needs of Alzheimer's patients. Through the integration of the PISTA Sound Tool, this approach leverages the science of brain entrainment to facilitate a self-questioning process aimed at achieving desired brain states.

The PISTA Sound Tool serves as a catalyst, stimulating the brain and inducing a state known as brain entrainment. This state is achieved when the brain synthesizes a third tone from two distinct tones delivered separately into the right and left ears. This alteration of sound frequencies renders the brain more amenable to change and constructive reinforcement. Notably, thought patterns undergo significant transformation. Unwanted and destructive thought patterns are systematically eliminated and replaced with constructive affirmations tailored to the specific cognitive and emotional challenges faced by Alzheimer's patients.

Within the framework of PISTA's holistic approach, seven distinct modalities are thoughtfully applied by a team of medical professionals with specialized expertise in psychology and mental health care. These professionals serve as guides, assisting patients in delineating their life direction and instilling the motivation to navigate the array of challenges encountered on a daily basis.

PISTA sessions are instrumental in helping both clients and therapists identify key triggers for adversities specific to Alzheimer's disease. These sessions facilitate the healing process by encouraging individuals to move forward in life through structured short and long-term programs, all under the compassionate and expert guidance of the dedicated PISTA team of professionals. The ultimate goal is to enhance mental well-being, promote resilience, and improve the overall quality of life for individuals grappling with the challenges of Alzheimer's disease.

PRINCIPLE OF NEUROPLASTICITY AND PERSONALIZED TRANSFORMATION

The methodology employed by PISTA is rooted in the principle of neuroplasticity and personalized transformation, with a focus on individuals living with Alzheimer's disease. It recognizes the brain's remarkable capacity to adapt, rewire, and optimize its functions even in the face of neurological challenges like Alzheimer's. This principle holds that by strategically leveraging a combination of brainwave feedback, brain stimulation, self-regulation techniques, and the power of sound entrainment through the PISTA Sound Tool, individuals can undergo profound cognitive and emotional transformations.

The methodology acknowledges that each individual's journey with Alzheimer's is unique, characterized by distinct cognitive, emotional, and physical challenges. Therefore, it emphasizes the importance of personalization, tailoring interventions to address the specific needs and goals of each patient. This personalized approach empowers individuals to embark on a journey of self-discovery, guided by a team of dedicated professionals with expertise in psychology and mental health care.

Through the principle of neuroplasticity and personalized transformation, PISTA strives to harness the brain's innate adaptability to replace destructive thought patterns with constructive affirmations, enhance cognitive function, and improve emotional well-being in individuals living with Alzheimer's disease. Ultimately, the methodology aims to optimize brain function, foster resilience, and improve the overall quality of life for those facing the challenges of Alzheimer's. It underscores the belief that, with the right interventions and support, positive and lasting change is not only possible but achievable for each individual on their unique path to cognitive and emotional well-being.

NEUROPLASTICITY

Neuroplasticity, often described as the brain's remarkable ability to reorganize itself, adapt, and form new neural connections, presents a profound opportunity for individuals living with Alzheimer's disease. This adaptive capability of the brain becomes particularly significant when considering the diverse challenges posed by Alzheimer's, which affect cognitive functions, memory, and emotional well-being. In Alzheimer's disease, the brain experiences the gradual loss of neurons and synaptic connections, resulting in cognitive decline, memory impairments, and emotional disturbances. Neuroplasticity intervenes by enabling the brain to compensate for these losses through various mechanisms. It allows the brain to create new neural connections, strengthen existing synapses, recruit healthy brain regions, and regulate emotions more effectively. By understanding and harnessing the brain's adaptive potential, individuals living with Alzheimer's can experience improved cognitive functions, better memory,

and enhanced emotional well-being. This approach provides a holistic and empowering perspective on managing the challenges associated with Alzheimer's disease.

PISTA'S INNOVATIVE APPROACH TO EMPOWER ALZHEIMER'S PATIENTS

PISTA, with its innovative approach tailored for individuals living with Alzheimer's, capitalizes on the principles of neuroplasticity and binaural beats to provide profound benefits. Neuroplasticity, the brain's remarkable capacity to adapt and rewire itself, takes center stage in enhancing the lives of Alzheimer's patients. PISTA harnesses the power of binaural beats to target working memory capacity by modulating brain activity within the 7.5 Hz – 12.5 Hz frequency range, recognized for its impact on a range of cognitive functions.

Extensive research supports the effectiveness of binaural beats in enhancing attention, inhibitory processes, perceptual abilities, and working memory. PISTA introduces the concept of the PISTA Wave, a carefully curated fusion of frequencies commencing at 7.5 Hz. These frequencies are thoughtfully selected to align with specific brainwave patterns, particularly those associated with cognitive functions like attention and working memory.

Additionally, the integration of PISTA Vibration adds an extra dimension to the program. Vibrations are employed in tandem with sound and speed to generate precise frequencies that harmonize with the program's objectives. This dynamic synergy of sound, speed, and vibration serves as a potent tool in facilitating the accomplishment of desired tasks and objectives, heightening the overall efficacy of the program.

By embracing these cutting-edge scientific principles and techniques, PISTA endeavors to unlock the full potential of neuroplasticity, extending the opportunity for cognitive enhancement, memory improvement, and enhanced cognitive functioning to individuals living with Alzheimer's. Through the strategic use of binaural beats, the PISTA Wave, and PISTA Vibration, this holistic approach represents a significant advancement in their journey toward an improved quality of life and cognitive well-being.

THE IMPACT OF PISTA AUDITORY BEAT STIMULATION

The PISTA approach, uniquely designed for individuals living with Alzheimer's, revolves around auditory beat stimulation to significantly impact their well-being. This comprehensive

strategy integrates various sound and vibrational models, crafting personalized modalities and protocols that seamlessly integrate into the daily lives and routines of Alzheimer's patients.

Auditory beat stimulation, a powerful cognitive tool, operates by delivering rhythmic auditory and vibrational patterns that generate perceptible "beats," actively engaging the brain. In the context of Alzheimer's disease, characterized by cognitive challenges and mood fluctuations, auditory beat stimulation emerges as a promising avenue for enhancing cognitive functions and mood regulation.

What sets PISTA apart is its commitment to empowering patients. Through the program, individuals with Alzheimer's do not merely receive therapy but actively participate in managing their conditions. PISTA equips patients with specific protocols during their coaching course, imparting valuable skills that enable them to navigate challenging situations and bolster self-confidence. Importantly, these skills extend beyond the coaching course, offering enduring support as patients confront and manage the cognitive and emotional complexities associated with Alzheimer's disease. PISTA represents a holistic, adaptable, and empowering approach to improving the lives of Alzheimer's patients by harnessing the potential of auditory beat stimulation.

PISTA'S EVIDENCE-BASED APPROACH FOR ENHANCING THE LIVES OF ALZHEIMER'S PATIENTS

The effectiveness of PISTA is substantiated by a wealth of research findings, including extensive longitudinal studies spanning three generations. These studies offer compelling evidence of the efficacy of PISTA's modalities, which encompass the application of sound and vibration across a spectrum of experiences, from memory and learning enhancement to addressing cognitive decline, managing panic attacks, and mitigating the effects of traumatic experiences. PISTA's unique approach within a psychological framework encourages the brain to process memories and emotions effectively, resulting in deep insights among research participants. Guided by trained PISTA coaches, individuals have harnessed these insights to instigate crucial changes, combat anxiety, reverse cognitive decline, and reframe traumatic experiences.

Crucially, the PISTA program integrates these research findings by harnessing the power of binaural beats and vibrations to establish a baseline of brain stimulation, particularly for working memory improvement. This approach exhibits considerable potential in benefiting individuals facing dementia and Alzheimer's disease, and it can also be effectively applied to

the daily lives of Parkinson's and Alzheimer's patients. Whether it's a leisurely walk in the park, a serene stroll by the sea, or the practice of meditation to cultivate a meditative state of mind, PISTA's integrated approach energizes and fosters motivated thinking for relevant tasks.

PISTA's rhythmic methodology strategically combines vibration speed and volume to achieve specific desired outcomes, inducing distinctive brainwave patterns that effectively alter negative neural connections. Through the integration of various frequencies, speeds, and volumes, PISTA actively regulates mood, enhances learning capabilities, bolsters memory, and supports cognitive and executive functioning. The program incorporates a variety of activities, including running, walking, and everyday tasks, to facilitate the entrainment of brainwave patterns.

What sets PISTA apart is its user-centric approach, allowing individuals to personalize programs based on their daily activities and desired effects. This model encompasses a wide range of rhythms, hertz frequencies, and machine models to cater to diverse needs, empowering individuals to embark on a personalized journey towards cognitive enhancement and improved overall well-being.

PISTA EMOTIONAL REGULATION

PISTA, meticulously tailored for individuals grappling with Alzheimer's disease, presents a multifaceted approach to regulating and improving emotional well-being. Alzheimer's often ushers in a host of emotions, including stress, fear, anger, despair, and a pervasive sense of demotivation. PISTA addresses these emotional challenges head-on, utilizing the therapeutic power of auditory beat stimulation to create a holistic impact.

One of the foremost emotional hurdles faced by Alzheimer's patients is stress. The unpredictability and cognitive decline associated with the disease can induce heightened stress levels. PISTA intervenes by incorporating rhythmic auditory and vibrational patterns, inducing a profound sense of calm and tranquility. This soothing effect aids in the management of stress, enabling patients to navigate the daily uncertainties with greater emotional resilience.

Fear and anxiety are common companions on the Alzheimer's journey, often stemming from the sense of confusion and loss of control. PISTA's approach offers a sanctuary of emotional stability. By fostering a serene environment through auditory beat stimulation, it imbues patients with feelings of security and comfort. This emotional steadiness empowers individuals to cope with the daunting uncertainties of Alzheimer's, providing a vital emotional buffer.

Anger and frustration can also manifest as emotional responses to the challenges posed by Alzheimer's. PISTA intervenes by promoting a sense of composure and balance. The rhythmic auditory cues and vibrations encourage a state of emotional equilibrium, helping patients approach adversities with a more level-headed perspective.

Alzheimer's can cast a shadow of hopelessness, leading to a profound lack of motivation. PISTA combats this by instilling hope and agency. The program equips patients with valuable skills, encouraging them to actively participate in their therapy. This sense of control and purpose can reignite the spark of optimism, bolstering their emotional well-being.

Moreover, PISTA enhances motivation by making therapy enjoyable and engaging. The auditory and vibrational elements transform the therapeutic process into an immersive experience, stimulating the mind and invigorating the spirit. This newfound motivation extends beyond therapy sessions, potentially improving overall engagement in daily activities and interactions.

In summary, PISTA's comprehensive approach caters to the intricate emotional landscape of Alzheimer's patients. By leveraging auditory beat stimulation, it addresses stress, fear, anger, hopelessness, and demotivation, fostering emotional resilience and enhancing overall quality of life for those navigating the complex terrain of Alzheimer's disease.

“WHY ME” Behavior

PISTA plays a pivotal role in addressing the "why me" behavior often exhibited by individuals living with Alzheimer's disease. This behavioral pattern, characterized by feelings of confusion, frustration, and questioning, can significantly impact the emotional well-being of patients. PISTA intervenes by providing a structured and supportive framework to help individuals better understand and cope with their condition.

One of the primary ways PISTA addresses this behavior is through cognitive restructuring. The program assists Alzheimer's patients in reframing their thoughts and perceptions. It empowers them to recognize that their condition is not a result of personal shortcomings but rather a complex neurological process. This cognitive shift helps individuals move away from self-blame and feelings of victimization, fostering a more accepting and self-compassionate mindset.

Additionally, PISTA leverages auditory beat stimulation to enhance cognitive functions such as memory and attention. By sharpening cognitive abilities, the program equips individuals to engage more effectively with their surroundings and daily routines. This heightened cognitive function can alleviate some of the distress associated with the "why me" behavior, as patients regain a sense of mastery and control.

Furthermore, PISTA offers emotional support and a sense of community. Group therapy sessions and interactions with trained PISTA coaches allow Alzheimer's patients to share their experiences and emotions with others who understand their struggles. This sense of belonging and empathy can mitigate feelings of isolation and helplessness, often associated with the "why me" behavior.

The personalized nature of PISTA is another key aspect of addressing this behavior. Tailored therapy plans and interventions consider each individual's unique challenges and needs. This individualized approach empowers Alzheimer's patients to take an active role in managing their condition, reducing feelings of powerlessness and frustration.

In summary, PISTA's holistic approach addresses the "why me" behavior in Alzheimer's patients by promoting cognitive restructuring, enhancing cognitive functions, providing emotional support and community, and tailoring interventions to individual needs. By fostering understanding, self-compassion, and empowerment, PISTA helps individuals navigate the emotional challenges of Alzheimer's disease with greater resilience and acceptance.

PISTA on Managing Fears and Anxiety

PISTA is a highly effective approach in managing fears and anxiety among individuals living with Alzheimer's disease. Alzheimer's can bring about significant emotional distress, including fear and anxiety stemming from confusion, memory loss, and uncertainty about the future. PISTA offers a multifaceted strategy to address these emotional challenges and enhance the well-being of Alzheimer's patients.

First and foremost, PISTA employs auditory beat stimulation to regulate mood and emotions. Through carefully designed binaural beats, the program can induce a calming effect on the brain, reducing stress and anxiety levels. These rhythmic auditory patterns have the potential to soothe the mind, providing relief from the persistent worries and fears commonly experienced by Alzheimer's patients.

Moreover, PISTA incorporates mindfulness and relaxation techniques as part of its therapy sessions. These practices help individuals develop a greater sense of emotional control and resilience. By teaching patients how to stay present in the moment and manage their responses to stressors, PISTA equips them with valuable tools for mitigating anxiety and fear.

Additionally, the personalized nature of PISTA interventions plays a crucial role in addressing fears and anxiety. Tailored therapy plans consider each patient's unique challenges and

emotional triggers. This individualized approach allows Alzheimer's patients to receive targeted support and coping strategies that address their specific concerns, ultimately reducing anxiety levels.

In conclusion, PISTA's holistic approach to managing fears and anxiety in Alzheimer's patients combines auditory beat stimulation, mindfulness practices, emotional support, and personalized interventions. By fostering a sense of calm, resilience, and community, PISTA empowers individuals to navigate the emotional challenges of Alzheimer's disease with greater confidence and emotional well-being.

PISTA on Depressive thoughts

PISTA, with its comprehensive approach, plays a pivotal role in facilitating a transformation in the patterns of depressive thought often experienced by individuals grappling with the challenges of Alzheimer's disease. This shift in thought patterns is achieved through a multi-faceted and holistic strategy that addresses the emotional and cognitive aspects of depression.

Firstly, PISTA incorporates cognitive restructuring techniques. Through cognitive training exercises, patients learn to identify and challenge negative thought patterns that contribute to depression. They acquire skills to reframe these thoughts, replacing them with more positive and realistic perspectives. This cognitive reprogramming is essential for breaking the cycle of depressive thinking.

Emotional processing is another cornerstone of PISTA's approach. Patients are provided with a safe and supportive environment to express and process their emotions, including the profound sadness and despair that often accompany Alzheimer's disease. This emotional expression is instrumental in acknowledging and addressing the emotional impact of the condition, allowing patients to work through their depressive thoughts.

Mindfulness and meditation practices are integrated into the PISTA program. These techniques encourage individuals to stay in the present moment, fostering a sense of inner calm and peace. They serve as effective tools for countering the persistent rumination and negative thought patterns that characterize depression.

Positive reinforcement is a key element of PISTA. Through objective progress tracking, patients receive regular feedback on their physical and cognitive improvements. This positive feedback serves as a powerful counterbalance to depressive thoughts, instilling a sense of accomplishment and hope for the future.

Furthermore, PISTA recognizes the importance of social support. Group therapy and support networks are encouraged, providing individuals with opportunities to share experiences and receive understanding and encouragement from peers who also face the challenges of Alzheimer's. This sense of community combats feelings of isolation and contributes to improved emotional well-being.

Building emotional resilience is a fundamental aspect of PISTA. Patients learn how to adapt to the realities of living with Alzheimer's disease, acquiring coping strategies and emotional fortitude. This resilience equips them with the tools to effectively manage depressive thoughts when they arise.

Ultimately, PISTA empowers individuals to take an active role in managing their condition. Through the practical tools and coping strategies, it provides, patients gain a sense of control over their health. This empowerment counters feelings of helplessness and reinforces the belief that they can influence their emotional well-being.

In conclusion, PISTA's approach to changing patterns of depressive thought is both comprehensive and holistic. By addressing cognitive restructuring, emotional processing, mindfulness, positive reinforcement, social support, resilience-building, and empowerment, PISTA empowers individuals to confront and overcome depressive thought patterns, fostering emotional well-being and a more positive outlook on life despite the challenges of Alzheimer's disease.

BUILDING CONFIDENCE WITH PISTA

PISTA plays a significant role in building confidence among patients with Alzheimer's disease, instilling in them a belief that their physical condition can improve. This confidence-building process is multifaceted and highly effective, thanks to PISTA's unique approach.

1. Objective Progress Tracking: PISTA employs objective measures to track the physical progress of patients. By regularly assessing and quantifying improvements in mobility, balance, and motor skills, patients receive tangible evidence of their physical advancements. This data-driven approach provides a sense of accomplishment and reinforces the belief that their condition can improve.

2. Physical Rehabilitation: PISTA integrates physical rehabilitation exercises tailored to the individual needs of each patient. These exercises are designed to enhance mobility, strength,

and coordination. As patients actively engage in these exercises and witness their own physical growth, their confidence in their ability to improve their condition naturally increases.

3. Cognitive Empowerment: Alzheimer's often involves cognitive challenges, and PISTA addresses this aspect as well. Through cognitive training exercises, patients sharpen their mental faculties, enhancing their problem-solving skills and cognitive flexibility. The resulting improvement in cognitive function bolsters their self-esteem and confidence in managing their condition.

4. Emotional Support: Confidence is closely tied to emotional well-being. PISTA creates a supportive environment where patients can openly express their fears and doubts. Through therapy sessions, they learn effective coping mechanisms and resilience-building techniques. This emotional support enhances their self-assurance in facing the physical challenges posed by Alzheimer's.

5. Empowering Autonomy: PISTA's approach is empowering, encouraging patients to take an active role in their own rehabilitation. By teaching patients how to manage their condition and providing them with tools and strategies, PISTA fosters a sense of autonomy and control over their health. This sense of agency significantly boosts their confidence.

6. Setting and Achieving Goals: PISTA encourages patients to set achievable goals for their physical and cognitive improvement. As they work towards and attain these milestones, their self-efficacy grows. The satisfaction of reaching these objectives reinforces the belief that further progress is possible.

In conclusion, PISTA's holistic approach to building confidence in Alzheimer's patients rests on a foundation of objective progress tracking, physical rehabilitation, cognitive empowerment, emotional support, autonomy, and goal-setting. Through these interwoven components, PISTA empowers patients to believe in their capacity for improvement, enhancing their confidence and overall quality of life.

PISTA Monitoring Progress and Personalized Support

PISTA devices play a pivotal role in monitoring the progress of Alzheimer's patients, providing a comprehensive assessment of how the physical illness impacts their daily lives. Through a multifaceted approach, these devices collect objective data on various aspects of the patient's physical and cognitive functioning, including gait patterns, motor skills, balance, and cognitive

performance. This data-driven approach allows healthcare providers to gain quantifiable insights into the patient's daily capabilities and limitations, offering a detailed view of the disease's impact. Additionally, patients can log their daily activities and experiences, creating a valuable record of their journey and highlighting any changes in symptoms or abilities. Longitudinal tracking over time enables clinicians to identify trends in the disease's progression, aiding in the adjustment of treatment plans. Regularly reviewed during in-person assessments with certified PISTA Coaches, this wealth of information aids in customizing interventions tailored to the patient's specific needs, optimizing the effectiveness of care and support for those living with Alzheimer's disease.

PISTA Devices

1. PISTA Power

Packed with a louder tone and stronger vibration, the PISTA Power Device is your ideal companion for supercharging your daily activities. Whether it's kickstarting your morning routine or maintaining peak productivity throughout the day, this device is designed to provide that extra push you need to elevate your performance and make the most out of every moment.



2. PISTA 24

The PISTA 24 Device is your essential companion for restful nights. Its soft tones and gentle vibrations create a serene ambiance, perfect for winding down your nighttime routine. Enjoy a tranquil transition to sleep and wake up refreshed, ready to conquer the day.



3. PISTA Sonic

PISTA Sonic is designed for seasoned PISTA users seeking advanced benefits and faster effectiveness. Elevate your wellness journey with these specially crafted headphones, enhancing your overall PISTA experience for accelerated results using sound and vibrations.



4. PISTA Life Power

This innovative device uses binaural beats and vibration technology to provide a truly immersive and therapeutic experience. With built-in vibration modes and heating pads that can be placed on the neck.



PISTA Beeps Speed and Frequency

The PISTA Device's binaural sound beats feature offers a range of different frequencies and speeds, allowing for a highly customizable experience tailored to an individual's specific needs. These binaural beats are essentially auditory stimuli created by playing two slightly different frequencies in each ear. When the brain processes these frequencies, it perceives a third frequency, known as the binaural beat, which is the difference between the two played frequencies.

The key feature of PISTA's binaural sound beats is their flexibility. Users can adjust the frequencies and speed of these beats to create a unique auditory experience. The frequencies can vary from slow to fast, offering a spectrum of options. This customization is essential because different frequencies and speeds can have distinct effects on the brain and cognition.

For example, slower beeps are associated with relaxation, stress reduction, and deep meditation. In contrast, faster beeps can enhance alertness, focus, and cognitive function. By allowing users to set the desired frequencies and speed, PISTA empowers individuals to choose the auditory stimulation that aligns with their specific goals and needs.

This customization is particularly valuable in the context of Alzheimer's disease, where cognitive and emotional well-being are crucial. Patients can use the PISTA Device to select the binaural sound beats that best suit their current state, whether it's to promote relaxation, boost cognitive function, or achieve a specific mental state. This adaptability ensures that individuals have access to a versatile tool for optimizing their mental and emotional health as they navigate the challenges of Alzheimer's.